

SOURCE CODE

TEAM ID	PNT2022TMID39429
PROJECT TITLE	PERSONAL ASSISTANCE FOR SENIORS WHO ARE SELF RELIANT

TINKERCAD CODE FOR MEDICAL REMAINDER

**// PRESS PUSH BUTTONS FOR TEMPERATURE AND
MEDICINE TIME REMAINDER**

#include <LiquidCrystal.h>

LiquidCrystal lcd(13, 12, 11, 10, 9, 8);

float voltage;

int celsius;

int valPulse;

int pulsePin = A2;

const int temperaturePin = A0;

String seconds;

long duration;

int buzzerpin = 2,button_1=A3,button_2 = 4;

void setup()

```
{  
  pinMode(temperaturePin, INPUT);  
  pinMode(button_1, INPUT);  
  pinMode(button_2, INPUT);  
  lcd.begin(16, 2);  
  Serial.begin(9600);  
  lcd.print("MEDICAL REMAINDER");  
  delay(2000);  
  lcd.clear();  
  lcd.setCursor(0, 0);  
  lcd.println("HEALTH MONITORING");  
  lcd.setCursor(0, 2);  
  lcd.print("SYSTEM");  
  delay(3000);  
  lcd.clear();  
  pinMode(buzzerpin, OUTPUT);  
}  
  
void loop()  
{  
  if(digitalRead(button_1) == HIGH)
```

```
{  
  lcd.clear();  
  lcd.print("Calculating....");  
  delay(3000);  
  voltage = analogRead(temperaturePin) * 0.004882814;  
  celsius = (voltage - 0.5) * 100.0;  
  lcd.clear();  
  lcd.print("Body Temp. : ");  
  lcd.print(celsius);  
  lcd.print("C");  
  delay(3000);  
  if(celsius>38)  
  {  lcd.clear();  
      digitalWrite(buzzerpin, HIGH);  
      lcd.print("High temp");  
      delay(2000);  
      lcd.clear();  
      lcd.setCursor(0, 0);  
      lcd.println("TAKE EMERGENCY");  
      lcd.setCursor(0, 2);
```

```
    lcd.print("TIME MEDICINE");  
    delay(3000);  
    lcd.clear();  
    digitalWrite(buzzerpin, LOW);  
}  
else  
{  
    lcd.clear();  
    lcd.print("Normal Temp");  
    delay(4000);  
}  
lcd.clear();  
lcd.print("Temperature:");  
lcd.print(celsius);  
lcd.print("C");  
lcd.setCursor(0,1);  
lcd.print("Pulse : ");  
lcd.print(valPulse);  
while(digitalRead(button_1)== 0);  
delay(250);
```

```
    lcd.clear();  
}  
if(digitalRead(button_2)== HIGH)  
{  
    while(1){  
        Serial.println("(hours) : ");  
        while (Serial.available()==0){}  
        long int hours = Serial.parseInt();  
        lcd.print(hours);  
        lcd.print(" : ");  
        Serial.println("(minutes) : ");  
        while (Serial.available()==0){}  
        long int mins = Serial.parseInt();  
        lcd.print(mins);  
        lcd.print(" : ");  
        Serial.println("(seconds) : ");  
        while (Serial.available()==0){}  
        long int seconds = Serial.parseInt();  
        lcd.print(seconds);  
        lcd.clear();  
    }
```

```
long int current_time =(hours*3600)+(mins*60)+(seconds);

    lcd.print(current_time);

    delay(3000);

    lcd.clear();

long int i;

for(i= current_time ; i< (current_time+(86400)); i++){

if(i == 28800){

digitalWrite(buzzerpin,HIGH);

tone(buzzerpin, 100 );

lcd.print("Medicine time");

delay(5000);

lcd.clear();

digitalWrite(buzzerpin,LOW);

}

else if(i == 72000){

digitalWrite(buzzerpin,HIGH);

tone(buzzerpin, 100);

lcd.print("Medicine time");

delay(5000);
```

```
lcd.clear();  
digitalWrite(buzzerpin,LOW);  
}  
else{  
    lcd.setCursor(0,0);  
  
    lcd.print("Your medicine");  
    lcd.setCursor(0,1);  
    lcd.print("time is after :");  
    delay(1000);  
    if(i<28800 ){  
        lcd.clear();  
        lcd.setCursor(0,0);  
        lcd.print(28800-i );  
        lcd.setCursor(0,1);  
        lcd.print("seconds");  
        delay(1000);  
        lcd.clear();  
    }  
    else if (i>28000 && i<72000){
```

```
    lcd.clear();  
    lcd.setCursor(0,0);  
    lcd.print( 72000-i);  
    lcd.setCursor(0,1);  
    lcd.print("seconds");  
    delay(1000);  
    lcd.clear();  
}  
else if (i>72000){  
    lcd.clear();  
    lcd.setCursor(0,0);  
    lcd.print( 86400-i + 28800);  
    lcd.setCursor(0,1);  
    lcd.print("seconds");  
    delay(2000);  
    lcd.clear();  
}  
}  
}  
}
```



```
}
```

```
}
```

PYTHON CODE FOR MEDICAL REMAINDER

```
import json
```

```
import wiotp.sdk.device
```

```
import time
```

```
import random
```

```
myConfig = {
```

```
    "identity": {
```

```
        "orgId": "tboyb4",
```

```
        "typeId": "medicineremainder",
```

```
        "deviceId": "19171603"
```

```
    },
```

```
    "auth": {
```

```
        "token": "12345678"
```

```
    }
```

```
}
```

```
client = wiotp.sdk.device.DeviceClient(config=myConfig,
logHandlers=None)

client.connect()

for i in range(0,20):

    tablet=["Amlodipine
Besylate","Azithromycin","Metformin","Amoxicillin","Cetiri
zine"]

    med icine time=[12.00,1.00,2.00,3.00,5.00,18.00,20.00,7.00]

    name = "Rajammal"

    medicine=random.choice(tablet)

    med icine time=random.choice(medicine time)

    mydata = {'Patient Name': name, 'Medicine Name':
medicine, 'Time': medicine time}

    client.publishEvent("IoT Sensor", "json", data=mydata,
qos=0, onPublish=None)

    print("Data published to IBM IOT platform :", mydata)

    time.sleep(5)

client.disconnect()
```